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Please read the usage instructions carefully before first use!

Icon legend

	Mark of conformity		max. weight load		Follow instruction for use
	Warning		Manufacturer		For single use only
	Keep dry		Keep away from sunlight		Marking Medical Device
	Item number		Batch / Lot no.		Date of production
	Quantity of handles		Fluid intake		Marking UDI

1. Technical data

REF	Article description	Dimension	Load capacity	Fluid intake	GTIN
4700181K	Multi- Umbettungstuch	80 x 210 cm	max. 200 kg	up to 4 Liter	Single pack (01)04050856047274 (11)XX (10)XX
					Shipping carton (01)04050856047281 (11)XX (10)XX

- 1.1. Material: Polyethylene film (blue), polypropylene fleece, superabsorbent, fluff pulp
- Material properties: Fluid resistance, Latex- and PVC-free, hypoallergenic/ dermatologically tested, x-ray compatible
- 1.2. Storage and transport
- Storage conditions: -5 bis 40°, protect from moisture and direct sunlight
- Humidity: k.A.

2. Indication / Areas of application / Intended use

The transfer sheet is an aid for transferring patients from one stretcher or bed to another (lifting, pushing or pulling). It also serves as an absorbent underlay for stretchers, stretchers and beds. The transfer sheet is primarily used in medical facilities, in the emergency services and in nursing care. It is not suitable for carrying patients over long distances. It is a disposable product.

3. Contraindications / risks of use

There are no known contraindications. If users are not trained or familiar with the correct use of the transfer sling, the sling may be used incorrectly. In addition, movements that are too fast and jerky can cause damage to the health of the person being moved and/or the user. It is imperative that four people are present to ensure safe transfer and temporary lifting. Observe the instructions for use under point 4 and the safety instructions under point 6. Open wounds or areas of skin that are not intact must be treated in accordance with the current recommendations for aseptic wound treatment and covered with a sterile dressing. It is the responsibility of the user to decide whether the dressing is suitable for the particular application. Transatlantic accepts no responsibility whatsoever for improper use of the reburial sheet.

4. Application / Placement

1. The entire length of the sheet (blue side down) is placed under the person to be transferred.
2. Each user (at least 4 persons) grasps two handles (one per hand = 8 handles)
3. Persons to be transferred must only be transferred using recommended/recognized and suitable safety equipment.

5. Service life / life cycle

The transfer pad is intended for single use.

6. Safety instructions / function control

WARNING indicates a possible hazard that could result in personal injury or property damage.

- It is the responsibility of the user to check, based on the circumstances (weight / state of health / type of injury of the person to be transferred), whether another means of transportation/rescue may be more suitable.
- Open wounds or areas of skin that are not intact must be treated in accordance with the current recommendations for aseptic wound treatment and covered with a sterile dressing.
- Make sure that a sufficient number of users are available to ensure the transfer of the body, depending on the situation.
- The maximum weight limit of 200 kg must be observed.
- Damaged transfer sheets (e.g. torn handles) may no longer be used.
- The transfer sheet is intended for single use only.

7. Service life / life cycle

The transfer sheet is a Class I medical device and complies with EU Regulation 2017/745 of the European Parliament and of the Council of 05.04.2017 for medical devices (MDR).

As a manufacturer and distributor of medical devices, we are certified in accordance with DIN EN ISO 13485. It goes without saying that our medical devices have undergone the relevant conformity assessment procedure and therefore bear the CE mark.

Example labeling:

8. Warranty

Transatlantic shall not be liable for any incidental or consequential injury or property damage.

9. Disposal information

There are no special disposal regulations for the Multi transfer pad. It must be disposed of in accordance with national and local regulations.

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