

#017

ANKLE STABILIZER FOR INVERSION AND EVERSION CONTROL



CHARACTERISTICS AND FUNCTIONALITIES

Helps recover as fast and better as possible:

- Adjustable 8-shaped bandage to provide progressive support and inversion/eversion fixation
- 3 adjustable hook & loop straps to allow a perfect closure

Comfortable:

- Breathable

User friendly:

- Back openable socks (easy to wear in case of painful and swelling ankle)

INDICATIONS

- Distortions
- Tendonitis
- Luxations



AMBIDEXTROUS

| SIZE | S | M | L | XL | XXL |
|------|-------|-------|-------|-------|-------|
| (CM) | 21-26 | 26-30 | 30-36 | 36-40 | 40-45 |

COLOURS ●



100% Breathable



#019

ANKLE BRACE WITH TWO RIGID SHELLS



CHARACTERISTICS AND FUNCTIONALITIES

Helps recover as fast and better as possible:

- External anatomical designed rigid shells (for medial and lateral immobilization) to provide ankle stability
- Adjustable (3 position) under foot band, to maximize fixing precision
- Inner air padding – for maximum support & comfort.

They also provide a micro massage during flexo and extension movements

- 2 adjustable closing bands – for adequate fixation

Comfortable: Breathable

Hygienic: Bacteriostatic Radyarn woven inner coating – antibacterial

User friendly:

- Completely openable (easy to wear in case of painful and swelling ankles). Adjustable base to fit different sizes

INDICATIONS

- Tibio / Tarsic distortions
- Luxations



Comfortable air padding and finished with antibacterial Radyarn fabric.



Adjustable for three different fitting sizes.

AMBIDEXTROUS

ONE SIZE

COLOURS ●

Hitech
TISSUE **rad yarn**
BACTERIOSTATIC bacteriostatic



100% Breathable



Prevents bad odour



Durable



Bacteriostatic

