

#674

ABDOMINAL BINDER ANATOMICAL

#675

ABDOMINAL BINDER UNISEX



CHARACTERISTICS AND FUNCTIONALITIES

Helps recover as fast and better as possible:

- Double and adjustable crisscrossing frontal strap - to allow a correct, progressive and comfortable closure regardless of patient shape
- Two back stays to give light support to the lumbar district

User friendly:

- Unisex

Comfortable:

- Breathable and elastic fabric

Hygienic:

- 100% soft cotton in contact with the skin
- Hypoallergenic fabric

• THE MODEL HIGH 18 CM CAN BE USED AS A CHEST UNIQUE BAND TOO

INDICATIONS

- Abdominal contention
- Post abdominal surgical op.



abdomen
rib cage

Model High 18 CM



AVAILABLE HEIGHT: 18 / 24 / 30 CM (#675) - 27 CM (#674)

SIZE	S	M	L	XL	XXL	XXXL
(CM)	70 - 82	82 - 94	94 - 105	105 - 117	117 - 130	130 - 150

COLOURS

**100%
COTTON**
ON YOUR SKIN



100%
Breathable



(only white)



#674-675