

#053

ACL-PCL AND MCL-LCL INSTABILITY BRACE



- The brace is studied to keep in the physiological position the knee joint during extension and flexion movements
- Two anatomical and totally openable fixings (one on the calf and the other on the thigh) makes the brace easy to wear. They are made by a special "3D Spacer woven" extremely breathable and that maximize comfort
- 4 elastic straps, closing with a "vise" movement" allows an easy and efficient brace fixation
- The two shaped stays (medial & lateral) give lateral knee stability
Are equipped with polycentric joints with variable rotation center, to perfectly reproduce the natural knee movement
- Both extension and flexion can be limited or a desired knee angle can be fixed
- Open on the knee (to facilitate medication and avoid stressing sensible areas)



INDICATIONS

- ACL-PCL fracture and/or reconstruction
- MCL-LCL fracture and/or reconstruction
- Support in case of light or medium gonarthrosis
- Support during active/passive rehab phases



RIGHT OR LEFT - HEIGHT: 40 CM

SIZE	XS	S	M	L	XL	XXL	* XXXL
(CM)*	34 - 38	38 - 43	43 - 49	49 - 56	56 - 64	64 - 73	73 - 83

* 15 cm above the knee - * Size XXXL available only upon request

COLOURS



3D Spacer

rad yarn
functional textile



100%
Breathable



Prevents
bad odour



Durable



Bacteriostatic

