



#579

STOD POSTURE TRAINER



- Our “Dynamic Orthopaedic Technical Support” **FITS EVERYONE**

- The shoulder straps are crossed on the back through a Hook & Loop closing system straps loop that can be adjusted according to the user's




- To warranty EFFICENCY regardless of patient physiognomy, two different bands closing positions are possible: A or B

- The abdominal band length can be cut with scissors
- To maximise COMFORT, two shoulder straps with different widths are provided


INDICAZIONI

- Postural disorders: incorrect postures, poor postural habits that cause lumbar and/or neck pain
- Posture misalignment syndrome
- Non-structural scoliosis
- Adjuvant in postural back pain

- **BREATHABLE** and **BACTERIOSTATIC** thanks to Pavis textile



HEIGHT: 40 CM



SIZE	S	M	L
Distance from shoulder to umbilicus (CM)	75 - 95	95 - 115	115 - 135



COLOURS ●



3D Spacer



100%
Breathable



Bacteriostatic

