

#021

# ADJUSTABLE ELASTIC KNEE BRACE WITH HOLE



### CHARACTERISTICS AND FUNCTIONALITIES

Helps recover as fast and better as possible:

- Frontal hole, with soft silicon ring, to provide patella stability
- Free popliteal area - to avoid compression in this critical area and warranty normal blood circulation
- 2 elastic straps to adjust compression according to leg shape and recovery necessity
- Two removable lightweight and flexible lateral stays to give further lateral stability to the knee

User friendly:

- Ambidextrous
- Completely openable

Comfortable:

- Breathable

Hygienic:

- 100% cotton in contact with the skin

### INDICATIONS

- Patellar sub-luxation
- Patellar instability
- Protective and contention action
- Joint effusions
- Meniscus traumas



AMBIDEXTROUS - HEIGHT: 29 CM

SIZE	SMALL	REGULAR	EXTRA
(CM)*	35 - 40	40 - 48	48 - 58

\* 10 cm Above The Knee

### COLOURS

**100% COTTON**  
ON YOUR SKIN



100%  
Breathable



#022

# ADJUSTABLE ELASTIC KNEE BRACE WITH HOLE AND POLYCENTRIC STAYS



### CHARACTERISTICS AND FUNCTIONALITIES

Helps recover as fast and better as possible:

- Frontal hole, with soft silicon ring, to provide patella stability
- Free popliteal area – to avoid compression in this critical area and warranty normal blood circulation
- 2 elastic straps to adjust compression according to leg shape and recovery necessity
- Two removable lightweight and flexible lateral stays to give further lateral stability to the knee
- Two removable polycentric shaped stays (one medial and one lateral) to give lateral stability to the knee and support MCL-LCL

User friendly:

- Ambidextrous
- Completely openable

Comfortable:

- Breathable

Hygienic:

- 100% cotton in contact with the skin

### INDICATIONS

- Patellar sub-luxation
- Patellar instability
- Protective and contention action
- Joint effusions
- Meniscus traumas

AMBIDEXTROUS - HEIGHT: 29 CM

SIZE	SMALL	REGULAR	EXTRA
(CM)*	35 - 40	40 - 48	48 - 58

\* 10 cm Above The Knee

### COLOURS

**100% COTTON**  
ON YOUR SKIN



100%  
Breathable



leg  
knee