

EXIT-POCKET® WP

Waterproof device
for catheter
exit-site protection



Description

Device for exit-site protection. Waterproof pocket closed on four sides provided with adhesive hole to be fixed to the skin.

Intended use

Containment of the catheter lumens during patient's hygiene and care procedures: suitable for protection of all exit-sites / catheters.

Indication

The device is completely waterproof: if well-positioned it protects the exit-site from getting in contact with water, preventing microbial contamination.

Available Models

CODE	HOLE	VERSION	INDICATION	PACKS/BOX
S2005H/WP	ø 5 cm	Closing	Bath / Shower	50
S2011H/WP	10 x 9 cm	Closing	Bath / Shower	
S2013HL/WP	13 x 11 cm	Closing	Bath / Shower	
S2013H/WP	11 x 13 cm	Closing	Bath / Shower	
S2007H/WP	5,5 x 7 cm	Closing	Bath / Shower	25
S2007HL/WP	7 x 5,5 cm	Closing	Bath / Shower	



S2005H/WP



S2011H/WP



S2013HL/WP



S2013H/WP



S2007H/WP



S2007HL/WP