

#576

CLAVICLE SUPPORT



CHARACTERISTICS AND FUNCTIONALITIES

Helps recover as fast and better as possible:


- Criss-cross back support system to pull back the shoulders
- Front Hook & Loop closing system to allow a desired closure

Comfortable: Elastic woven, Soft padding under the armpit, Breathable

INDICATIONS

- Hypercyphosis
- Correction in case of shoulders anteposition

Detail of the crossing of the two bands on the dorsal part.



SIZE	S	M	L	XL	XXL	XXXL
(CM)*	70-76	76-82	82-88	88-96	96-102	102-110

* Distance from shoulder to umbilicus

COLOURS ○



100%
Breathable



#577

CLAVICLE IMMOBILIZER



CHARACTERISTICS AND FUNCTIONALITIES

Helps recover as fast and better as possible:


- Criss-cross back support system to pull back the shoulders
- Hook & Loop closing system straps

Comfortable:

- Soft padding under the armpit
- Breathable

INDICATIONS

- Clavicle fractures



SIZE	S	M	L
(CM)	65 - 85	85 - 105	105 - 125

COLOURS ●



100%
Breathable

